




QUICK START GUIDE

FINGERTIP PULSE OXIMETER

HOLFENRY

1



INSERT FINGER

Nail facing up,
insert completely



2



PRESS POWER

Press the button
to turn on the device



3



STAY STILL

Keep hand and body
still during measurement



4

READ RESULTS

View your SpO₂, Pulse Rate
and PI%

Blood Oxygen
Saturation(SpO₂)



Pulse Rate(PR) /
Perfusion Index(PI)

5

AUTO SHUTOFF

Remove your finger,
8 seconds auto
power-off



DISPLAY REFERENCE

Blood Oxygen
Saturation(SpO₂)



Pulse Rate(PR) /
Perfusion Index(PI)



Click the button
can show the
Pulse Rate or PI



Note: The fingernail should be facing the top chamber with sensor.
Please insert your finger completely into the chamber.



For best results, keep finger warm, clean and dry. Do not wear nail polish.



Only read data when the reading stays stable.

